

# Overhead Stirrer

---



## IKA RW 20 Digital Overhead Stirrer

The bestseller in the laboratory:

- With digital display
- Robust, slimline, ergonomic design
- Technical improvements on the trusted RW 20 series designs
- With constant power-drive
- Two speed ranges for universal use from 60 - 2,000 rpm
- Two speed ranges for universal use from 60 - 2,000 rpm
- Push-through agitator shafts (only when stationary)
- For stirring quantities of up to 20 L (H<sub>2</sub>O)

## Technical Data:

Stirring quantity max. per stirring position (H2O) [L]	20
Speed range [rpm]	60 - 2000
Viscosity max. [mPas]	10000
Speed control	Stepless
Dimensions (W x H x D) [mm]	88 x 294 x 212
Weight [kg]	3.1
Permissible ambient temperature [°C]	5 - 40

[www.severnsaleslabequip.com](http://www.severnsaleslabequip.com)



+44 (0) 117 9608858  
[sales@severnsaleslabequip.com](mailto:sales@severnsaleslabequip.com)