

B-ONE PLUS VISCOMETERS



Intuitive and powerful, the new B-ONE PLUS will surprise you with its extended measuring range in LR version or with its optional temperature sensor.

- 7" Touch Screen.
- Easy attachment of your spindle.
- Available in standard and high sensitivity versions.

Rack stand in option



Portable B-ONE
with carry case



Specification:

Type of instrument	Rotating springless viscometer with 7" Touch screen
Rotation speeds	Unlimited number of speeds between 0.3 and 250 rpm
Torque range	Standard Version : 0.05 to 13 mNm LR Version : 0.005 to 0.8 mNm
Accuracy	+/- 1 % of the full scale
Repeatability	+/- 0,2 %
Display	Viscosity – Speed – Torque – Time - (Temperature in option) Choice of viscosity units: cP/Poises or mPa.s / Pa.s
Security and confidentiality	An «operator» function allows you to enter a username for your instrument. There is also a protected mode that locks your measurement conditions.
Language	French/English/Russian/Spanish
Compatible measuring system	MS ASTM, MS BV, MS VANE
Compatible temperature control	EVA LR-BV
Supply voltage	90-240 VAC 50/60 Hz
Options	Carry case (PN 200500) PT100 probe (-50oC to +300oC) (PN 900026) Rack stand (PN P008000)
Dimensions and weight	Head: L180 x W135 x H250 mm Hardened steel stand: L280 x W200 x H30 mm Stainless steel rod: Length 500 mm Weight: 6.7 kg

www.severnsaleslabequip.com



+44 (0) 117 9608858
sales@severnsaleslabequip.com